

## Daily Checklist for labour ward resuscitation room Resuscitaire #1

Daily clean and check all equipment and drugs including expiry dates. Clean and check after each use. Student midwives may check resuscitaires without supervision once deemed competent by mentor. Please ensure all stock with specified levels are correct and replace as necessary.

<b><u>Date:</u></b>								
<b><u>Resuscitaire</u></b>								
Air cylinder x 1								
Oxygen cylinder x 1								
Clock								
Check Oxygen hoses for wear and tear								
Check piped gases functioning following switch from cylinders back to piped gases								
Towels and linen								
Heater and light functioning								
<b><u>Ready for use on top of Resuscitaire</u></b>								
Laerdal bag with medium mask [01]								
ET Catheter Mount								
Suction tubing with medium connection & 10fg catheter or yankauer								
Stethoscope								
<b><u>Drawer under Resuscitaire – Top Tray</u></b>								
Infant Laryngoscopes – with blades + extra small blade								
ET Introducers x 2								
Stitch or disposable sterile scissors x 1								
Spencer Wells x 1								
Infant guedel airways 0 x 1 00 x 1 000 x 1								

Medium Connectors x 4								
NG tubes 6fg, 8fg x 1 each								

<b><u>Date:</u></b>								
<b><u>Drawer under Resuscitaire</u></b>								
Extension set with 'T' adaptor x 4								
Jelco 24G [yellow] IV cannula x 4								
Bandage								
Splints x 2								
Transpore 2.5 mm x 1								
Yale Spinal needles 22G x 3								
Green and Orange Needles								
Syringes 1 ml x 4, 2ml x 4, 5ml x 2, 10 ml x 2								
Pink Lancets – Premie								
50 ml Syringes & White needles x 2								
Tempadots								
Cord Clamps x 4								
PH Indicator Strips								
Naloxone 400mcg x 2 amps								
Sodium Chloride 0.9% 10 mls x 10								
Dextrose 10% x 1 500ml bag								
Adrenaline 1:10,000 syringe x 2								
Sterets								
Portex tracheal tube with holder and clamp 1 x each 2.0mm,2.5mm, 3.0mm, 3.5mm, 4.0mm								
2 x 10ml enteral feeding syringe								
<b><u>Resuscitation Trolley</u></b>								
<b><u>Drawer 1</u></b>								
Portex tracheal tube with holder and clamp 1 x each 2.0mm,2.5mm, 3.0mm, 3.5mm, 4.0mm								
Laryngoscope spare blades								
Infant guedel airways sizes 0, 00, 000 x 1 each								
Sterile Disposable scissors								
Wells Forceps								
<b><u>Drawer 2</u></b>								
Suction Catheters 6, 8, 10 fg								
Yankauer Suckers x 2								
Oral/enteral Feeding syringes 2ml, 5ml, 20ml x2 each								
PH Indicator Strips								
NG tubes size 6 and 8								



